

Application Guide:
Light Collimation Solutions with Mighty Light and FWS Poly


Technical details for collimated output from our Mighty Light broadband light source and FWS Poly tunable filter device.

[View Details](#)

Mighty Light & Collimator combination



Mighty Light + ML-CA-S + Poly(CW: 660 nm, FWHM: 10 nm)



Exit diameter	10 mm
Beam diameter (200 mm propagated)	25 mm
Beam divergence	2.14 °
Output power	103 μW

[See All Products](#)

Wavelength Selector

A truly revolutionary tool for spectroscopy and spectral imaging.

[View Details](#)



Light Sources

Spectrolight offers a growing selection of unique light sources and related accessories


[View Details](#)



Optical Components

Spectrolight provides a growing range of optomechanical components

[View Details](#)



Systems

Spectrolight provides a growing range of innovative systems and related accessories

[View Details](#)



For Sales and Marketing Inquiries

sales@spectrolightinc.com
949-800-7780

For Applications and Technical Inquiries

info@spectrolightinc.com

SLI
www.spectrolightinc.com



Copyright © 2018 Spectrolight Inc., All rights reserved.

Our mailing address is:
info@spectrolightinc.com

Want to change how you receive these emails?

You can [update your preferences](#) or [unsubscribe from this list](#).