

**LaserFocusWorld 2020 Innovator's  
Awards GOLD Honoree!**

We are very pleased to announce our Flexible Wavelength Selector for Hyperspectral Imaging (FWS-HSI) is a Gold honoree in Detectors & Imaging.




[View Details](#)

**Optical ruler for laser beam alignment**

The Light Aligner is an optical ruler with precision pinholes and fiducials that simplifies the common problem of aligning laser beams. A new bolt-on accessory for Infrared(IR) detection is also available.

[View Details](#)





As you can see, when IR slider comes into the path of IR laser beam it can be visualized as an orange spot.

### **Light Aligner IR Slider**

The IR slider is a bolt-on accessory of Spectrolight's beam alignment tool, Light Aligner. It can be used to visualize and align infrared beams. As you can see, when the IR slider comes into the path of an IR laser beam it can be visualized as an orange spot. The wavelength range for detection is from 700 to 1600 nm with a minimum intensity of around  $3\mu\text{W}/\text{cm}^2$ . The screw can be tightened to fix the location of the IR slider.

[View Details](#)

[See All Products](#)

### **Wavelength Selector**

A truly revolutionary tool for spectroscopy and spectral imaging.

[View Details](#)



### **Light Sources**

Spectrolight offers a growing selection of unique light sources and related accessories


[View Details](#)



## Optical Components

Spectrolight provides a growing range of optomechanical components

[View Details](#)



## Systems

Spectrolight provides a growing range of innovative systems and related accessories

[View Details](#)



For Sales and Marketing Inquiries

[sales@spectrolightinc.com](mailto:sales@spectrolightinc.com)  
949-800-7780

For Applications and Technical Inquiries

***Jun Hee Kang, Ph. D.***  
Application Scientist

Email: [jkang@spectrolightinc.com](mailto:jkang@spectrolightinc.com)  
Mobile: 949-800-5117



Copyright © 2018 Spectrolight Inc., All rights reserved.

**Our mailing address is:**  
[info@spectrolightinc.com](mailto:info@spectrolightinc.com)

Want to change how you receive these emails?  
You can [update your preferences](#) or [unsubscribe from this list](#).